

Quick Installation and Maintenance Instructions

Patient Monitor Mounting




MON-VESA-ADA-DRP-S

Medical Facility Responsibility


Preventive maintenance checks must be performed regularly to maintain the quality and performance of this product. Any parts that may be broken, missing, worn, distorted or contaminated in any way should not be used and all affected parts should be replaced immediately. Should the necessity of any repair be suspected; please contact your local Amico distributor.

⚠ WARNING: It is the responsibility of the end user to ensure all aspects of installation are covered.

Installation Tools



Installation Reference




Installation to Device (Spacelabs Qube/Qube Mini)

ⓘ NOTE: Insert and tighten screws in a cross pattern.

INSTRUCTIONS:


1. Lay monitor on a padded surface and align parts as shown.
2. Insert 4 x H-PHMS-M4-10 through X-MON-VESA-ADA-DRP-S and into the monitor.



Installation to Mounting Platform (Amico #2 Head and -S Head)

1. Pull out X-MON-PLUNPIN from X-DROP-SHLF and rotate 90 degrees to lock in outer position*.
2. Ensure 2 x H-PHSCS-1420-12 are loosened and do not protrude from the front surface of X-DROP-SHLF**.

*PULL OUT AND ROTATE PLUNGER PIN [X-MON-PLUNPIN] 90°



**NO PROTRUSION ON FRONT SURFACE OF SHELF [X-DROP-SHLF]


Quick Installation and Maintenance Instructions

Patient Monitor Mounting



MON-VESA-ADA-DRP-S

- With monitor mounted on X-MON-VESA-ADA-DRP-S (not shown), slide X-MON-VESA-ADA-DRP-S into X-DROP-SHLF by lining up 2 x X-MON-BLK-STDOF with the 2 x cutouts on X-DROP-SHLF* so that 2 x X-MON-BLK-STDOF and X-MON-VESA-ADA-DRP-S sandwich the thinner portion of the 2 x cutout on X-DROP-SHLF**.
- Rotate X-MON-PLUNPIN 90 degrees to release it through the corresponding hole in X-MON-VESA-ADA-DRP-S***, and hand tighten 2 x H-PHSCS-1420-12 to secure X-MON-VESA-ADA-DRP-S in X-DROP-SHLF.



Troubleshooting Guide

Symptom	Possible Cause	Solution
Monitor is not secure	Loose mounting to bracket	Re-align parts and refer to instructions (Installation to Device).
Monitor is secure but bracket is not	Loose mounting to mounting platform	Re-align parts and refer to instructions (Installation to Mounting Platform).

Maintenance

It is recommended that once every three (3) months, all locks, bolts and screws are visually inspected for signs of product wear and tear or damage. Check the product for looseness in the mount or orientation.

Warranty

During the term of warranty: Within the first twelve (12) months from the date of shipment, Amico Accessories will repair or replace any part which is proven to be defective at no cost. After the twelve (12) month period, Amico Accessories will send the parts to the customer free of charge, however shipping and installation will be borne by the customer.

The warranty is valid only when the product has been properly installed according to Amico Accessories specifications, used in a normal manner, and serviced according to factory recommendations. It does not cover failures due to damage which occur in shipments or failures which resulted from accidents, misuse, abuse, neglect, mishandling, alteration, misapplication or damage that may be attributable to Force Majeure.

AMICO ACCESSORIES DOES NOT HONOR VERBAL STATEMENTS CONCERNING THE WARRANTY.

The distributor and/or dealer are not sanctioned to create verbal warranties about the product described in this agreement. Any statements will not be honored or be made part of the agreement of sale. This document is the final complete and exclusive terms of the agreement.

THIS WARRANTY IS INCLUSIVE AND REPLACES ALL OTHER WARRANTIES.

Amico Accessories shall not, under any circumstances be liable for incidental or consequential damages including, but not limited to, profit, loss of sales or injuries to person(s) or property.

Correction of non-compliance as noted above will result in completion of all liabilities of Amico Accessories whether based on agreement, neglect or changed materials, designs or specifications without notice.

All claims for warranty must first be approved by Amico Accessories Customer Service Department: (info@amico-accessories.com or 1-877-264-2697). A valid Return Goods Authorization number must be obtained from Amico Accessories prior to commencement of any warranty claim.

Amico Accessories Inc. 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada
Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587

www.amico.com