








# Nihon Kohden Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Life Scope G9 LCD Display Unit 12.12 lbs (5.50 kg)*</p>	N/A	 <p>AHM-057</p>	N/A	 <p>Core Unit CSM-1900</p>  <p>BSM1700 + JA694PA</p>
		 <p>SSM-057</p>	N/A	 <p>WS-930P</p>  <p>Cable and Patient Leads</p>

\*Weight provided to help approximate mounting capacity. For reference only.

\*\*Vertical rail and anesthesia/mobile cart options available. Please inquire.

# Nihon Kohden Ordering Guide Reference Chart









Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p><b>Core Unit CSM-1900</b> 27.55 lbs (12.50 kg)*</p>	 <p>NK-CSM19-ADA-M</p>	N/A	Mounts on MRS rail with AHM/SSM	N/A
 <p><b>BSM1700 + JA694PA</b> BSM1700: 3.46 lbs (1.57 kg)* JA694PA: 4.40 lbs (2.00 kg)*</p>	 <p>NK-JA694PA-ADA-R125</p>	 <p>AHM-057</p>	Mounts below with primary device on AHM/SSM	N/A
	 <p>DOWNPOST-16</p>	 <p>SSM-057</p>		
 <p><b>WS-930P</b> 1.7 lbs (0.78 kg)*</p>	 <p>NK-QI930P-ADA-M</p>	N/A	Mounts on MRS rail with AHM/SSM	N/A

\*Weight provided to help approximate mounting capacity. For reference only.

\*\*Vertical rail and anesthesia/mobile cart options available. Please inquire.

# Nihon Kohden

## Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons	
 <p>Cable and Patient Leads</p>	 <p>CRDHLDR</p>	 <p>AHM-057</p>	Mounts below with primary device on AHM/SSM, may require a downpost	N/A	
		 <p>SSM-057</p>			
	 <p>CRDHLDR2</p>	 <p>AHM-057</p>	Mounts below with primary device on AHM/SSM, may require a downpost		
		 <p>SSM-057</p>			
	 <p>101004-MON-BB</p>	N/A	Mounts on MRS rail with AHM/SSM		N/A

\*Weight provided to help approximate mounting capacity. For reference only.

\*\*Vertical rail and anesthesia/mobile cart options available. Please inquire.