
















Spacelabs Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p data-bbox="275 976 491 1049">Spacelabs Qube 9 lbs (4.1 kgs)*</p>	 <p data-bbox="695 922 947 948">MON-VESA-ADA-DRP-5</p>	 <p data-bbox="1100 443 1203 469">AHM-05S</p>	N/A	 <p data-bbox="1759 591 1877 617">Capno Pod</p>
		 <p data-bbox="1104 662 1199 688">SSM-05S</p>	N/A	 <p data-bbox="1745 878 1885 904">Patient Leads</p>
		 <p data-bbox="1062 967 1241 993">TTM-05-45S-A00</p>	N/A	 <p data-bbox="1745 1227 1885 1253">Power Supply</p>
		 <p data-bbox="1083 1336 1224 1362">RS6-45S-PKG</p>	N/A	






*Weight provided to help approximate mounting capacity. For reference only.

Spacelabs Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Spacelabs Qube with Capno Pod Qube - 9 lbs (4.1 kg)* Capno Pod - 0.49 lbs (0.22 kg)*</p>	 <p>MON-VESA-ADA-DRP-5</p>	 <p>AHM-05S</p>	<p>(Mount to the back of the Spacelabs Qube)</p>	<p>N/A</p>
		 <p>SSM-05S</p>		
		 <p>TTM-05-45S</p>		
		 <p>RS6-45S-PKG</p>		

*Weight provided to help approximate mounting capacity. For reference only.

Spacelabs Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Patient Leads</p>	N/A	 <p>101004-MON-BB</p>	Mounts on MRS rail with AHM/SSM	N/A
 <p>Power Supply 1.5 lbs (0.90 kgs)*</p>	N/A	 <p>PWRHLD-MON</p>	Mounts on MRS rail with AHM/SSM	N/A
		 <p>PWRHLD-R125</p>	Use R125 adapter when mounting on roll stand or TTM	

*Weight provided to help approximate mounting capacity. For reference only.