















Spacelabs Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p data-bbox="249 976 514 1045">Spacelabs Xprezzon 9.1 lbs (4.3 kgs)*</p>	N/A	 <p data-bbox="1052 586 1266 610">AHM-LCD-05-01-SL</p>	N/A	 <p data-bbox="1745 467 1892 492">Xprezzon CPU</p>  <p data-bbox="1749 721 1885 745">Patient Leads</p>
	N/A	 <p data-bbox="1056 1141 1262 1166">SSM-LCD-05-01-SL</p>	N/A	 <p data-bbox="1745 1040 1892 1065">Power Supply</p>  <p data-bbox="1766 1325 1871 1349">XP-SL-XPZ</p>





*Weight provided to help approximate mounting capacity. For reference only.

Spacelabs Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Xprezzon CPU 9.1 lbs (4.3 kgs)*</p>	 <p>X-SL-XP293-LCD-CPUMNT-BR</p>	 <p>AHM-LCD-05-01-SL</p>	<p>Mount at the back of the Monitor</p>	<p>N/A</p>
		 <p>SSM-LCD-05-01-SL</p>		
	 <p>X-SL-XP293-LCD-CPUMNT-BR + DOWNPOST-03</p>	 <p>AHM-LCD-05-01-SL</p>	<p>Mount below with post primary device on AHM/SSM</p>	
		 <p>SSM-LCD-05-01-SL</p>		

*Weight provided to help approximate mounting capacity. For reference only.

Spacelabs Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Patient Leads</p>	 <p>XP-SL-XPZ</p>	N/A	Mount at the back of the Monitor	N/A
 <p>Power Supply 1.5 lbs (0.90 kgs)*</p>	N/A	 <p>PWRHLD-MON</p>	N/A	N/A

*Weight provided to help approximate mounting capacity. For reference only.