



Key Features:

- 550 lbs (250 kg) weight capacity
- 23" (58 cm) low height
- Storage area for patient belongings
- Radiolucent fowler section
- Easily accessible brake pedals
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces
- Single-hand-operated lowering leg section
- Retractable fluid basin holder

Standard Features:

- Two IV pole holders
- Roller bumpers
- Adjustable fowler section
- Fifth wheel steering
- Two side mounted pedals for patient surface height adjustment
- Hydraulic lift system capable of supporting up to 550 lbs (250 kg)
- Central brake locking system
- Ergonomic push handles
- Side bumpers
- Head and foot-end mattress retainer
- Retractable fluid basin holder
- Adjustable calf supports

The Amico Titan OB/GYN Stretcher is designed with functionality and convenience in mind, allowing caregivers to provide a more meaningful level of care to patients.

This product offers key features such as a one hand operated lowering leg section, a deep well storage area for patient belongings and an ergonomic design to help improve efficiency of healthcare. To accommodate obstetrics patients with limited mobility, the Titan OB/GYN Stretcher is built on Amico's best-in-class, low height stretcher platform.

The Titan OB/GYN Stretcher can also serve as an interim transport stretcher for moving patients between departments.

Technical Specifications

Model Number & Approval

Model Number	S-AM1-3005
Agency Approval	cULus

Overall Dimensions

Overall Length	84" (213 cm)
Overall Width	33" (84 cm)

Litter Positioning

Trendelenburg/Reverse Trendelenburg	16°
Maximum Fowler Section Inclination	80°

Height Range

Lowest	23" (58 cm)
Highest	33" (83.8 cm)

Surface Dimensions

Surface Width	25" (63.5 cm)
Surface Length	75" (191 cm)
Surface Thickness	4" (10 cm)

Weight Capacity & Caster Diameter

Weight Capacity	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- 4"/5" (10 cm/12.7 cm) high density foam mattress with or without fire barrier
- Disposable fluid basins
- Permanent IV pole
- Foley bag holder
- Overbed tables
- Stirrups (calf supports)

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: Three years parts