

Titan Series

Titan Transport



Key Features:

- 550 lbs (250 kg) weight capacity
- Standard fifth wheel for improved mobility
- 21" (53 cm) low height
- Storage area in base for two oxygen E-cylinder tanks and belongings
- Radiolucent fowler section
- Dual brake pedal mechanisms
- Optional knee gatch for enhanced comfort
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces

The Amico Titan Transport Stretcher uses an innovative design and incorporates distinctive features aimed to help caregivers provide a more meaningful level of care to patients.

This product is specifically designed to meet the needs of customers by offering key characteristics such as a dual brake system, deep well storage area in the base and an ergonomic design to help improve efficiency of healthcare. With the heightened awareness regarding fall prevention, Amico is proud to offer a Transport Stretcher that is the best-in-class low height stretcher in the market.

Standard Features:

- Four IV pole holders
- Roller bumpers
- Adjustable fowler section
- Fifth wheel steering
- Stationary foot end litter
- Fold-down siderails
- 8" (20 cm) wheel casters
- Hydraulic lift system
- Central brake locking system
- Ergonomic push handles
- Side bumpers
- Head and foot-end mattress retainer
- Integrated foot push handle

Model Number & Approval

Model Number	S-AM1-3000
Agency Approval	cULus

Overall Dimensions

Overall Length	84.5" (214 cm)
Overall Width	33" (84 cm)

Litter Positioning

Trend/Reverse Trend	17°
Maximum Fowler Section Inclination	80°
Maximum Knee Section Inclination	35°

Height Range

Lowest	21" (53 cm)
Highest	33" (84 cm)

Surface Dimensions

Surface Width	25" (63.5 cm)
Surface Length	75" (191 cm)
Surface Thickness	4" (10 cm)

Weight Capacity & Caster Dimension

Maximum Weight Limit	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- 4"/5" (10 cm/12.7 cm) high density foam mattress with or without fire barrier
- Overbed tables
- Foley bag holder
- Permanent IV pole
- Utility tray

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: Three years parts

Amico Beds Corporation 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 71 East Industry Court, Deer Park, NY 11729, USA
Toll Free Tel: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862