



Key Features:

- Easily removable head and foot boards
- Instant CPR release lever
- 550 lb (250 kg) weight capacity
- LED night light
- Bed bumpers for added protection
- Five positional functions (HILO, Fowler, Knee, Trendelenburg and Reverse Trendelenburg)
- Battery back-up
- Caregiver-controlled "flat surface" function
- Fowler angle indicator
- 6 IV pole holders
- 2 Foley bag holders
- Linak DC actuators
- Patient control on inner siderails

The Amico Apollo Hospital Bed Collection provides caregivers the control they need in the recovery environment. The simple design and distinctive features bring efficiency to today's medical and surgical environment.

The Amico Apollo AC 400C features easily removable head and foot boards which provide caregivers easy access in the event of an emergency. The head and foot boards feature an exclusive blend of polypropylene for uncomplicated cleanup and durability.

The functions of the bed are electronically controlled and equipped with a back-up battery that allows caregivers to make adjustments to the bed during transport and unexpected power outages.

Model Number

Model Number	AC-400C
--------------	---------

Overall Dimensions

Overall Length	88" (224 cm)
Overall Width	42" (107 cm)
Patient Surface	35" x 77" (89 x 196 cm)
Height Range (to the surface of the litter)	
High	29" (74 cm)
Low	18" (46 cm)

Litter Positioning

Trend./Reverse Trend.	$\pm 12^\circ$
Backrest	0° - 65°
Knee Gatch	0° - 30°

Weight Capacity & Caster Diameter

Weight Capacity	550 lbs (250 kg)
Caster Diameter	5" (13 cm)

Electrical

Voltage	110 V @ 60 Hz or 220 V @ 50 Hz
---------	--------------------------------

Optional Features & Accessories:

- Chart holder
- Provision for Trapeze
- 6" (15 cm) caster

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: Three years parts