














GE Healthcare Ordering Guide Reference Chart







Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Transport Pro Monitor 8.4 lbs (3.8 kg)*</p>	 <p>GE-TRPRO-BKT-00 (Allows GE PDM connection to the back of the monitor) 0.5 lbs (0.22 kg)*</p>	 <p>AHM-057</p>	N/A	 <p>RAC2</p>  <p>Power Supply</p>  <p>Cuffs and Patient Leads</p>  <p>GE Unity</p>  <p>GE PDM</p>
		 <p>SSM-057</p>	N/A	
		 <p>TTM-16-457-A00</p>	N/A	
	<p>N/A (The Transport Pro has a 100x100 VESA pattern on the back, but the RAC2 and GE PDM will not be able to mount to the back of the Transport Pro)</p>	 <p>AHM-057</p>	N/A	
		 <p>SSM-057</p>	N/A	
		 <p>TTM-16-457-A00</p>	N/A	

*Weight provided to help approximate mounting capacity. For reference only.

GE Healthcare Ordering Guide Reference Chart








Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>RAC2 5 lbs (2.3 kg)*</p>	N/A	 GE-RAC2M-MON-00	Mounts on MRS rail with AHM/SSM	N/A
	 GE-RAC4-SI-00 1.5 lbs (0.68 kg)*	 AHM-057  SSM-057	Mounts below with primary device on AHM/SSM	
 <p>Power Supply 1.5 lbs (0.90 kg)*</p>	N/A	 PWRHLD-MON	N/A	N/A
 <p>Cuffs and Patient Leads</p>	N/A	 100609-MON-BB	Mounts on MRS rail with AHM/SSM	N/A

*Weight provided to help approximate mounting capacity. For reference only.

GE Healthcare Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 GE Unity 2.4 lbs (1.1 kg)*	 GE-UNITY-PLT-00 0.28 lbs (0.12 kg)*	 GE-8000-MON-00	N/A	N/A
 GE PDM 2.8 lbs (1.3 kg)*	N/A (The GE PDM mounts to the back of the Transport Pro only with GE-TRPRO-PLT-00)	N/A	N/A	N/A
	N/A	 GE-PRNPD-LBT-MON	Should not be used with GE-TRPRO-BLT-00	N/A

*Weight provided to help approximate mounting capacity. For reference only.