









Philips Ordering Guide Reference Chart





Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Philips IntelliVue MP20/30 13.4 lbs (6 kg)*</p>	 <p>N/A (Philips supplied Tabletop Mounts)</p>	 <p>AHM-053</p>	N/A	 <p>Multi-Measurement Module (MMS)</p>  <p>Cable and Patient Leads</p>
		 <p>SSM-053</p>	N/A	
		 <p>TTM-05-453-A00</p>	N/A	
		 <p>RS6-423-PKG</p>	N/A	

*Weight provided to help approximate mounting capacity. For reference only.

Philips Ordering Guide Reference Chart









Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Philips IntelliVue MP20/30 13.4 lbs (6 kg)*</p>	 <p>N/A (Philips supplied M8045A Docking Station with Tabletop Mount has VESA 100x100) 5.95 lbs (2.7 kg)*</p>	 <p>AHM-050</p>	N/A	 <p>Multi-Measurement Module (MMS)</p>  <p>Cable and Patient Leads</p>
	 <p>MON-VESA-ADA</p>	 <p>SSM-050</p>	N/A	
	 <p>N/A (Philips supplied M1026A/B, M1013A or M1019A Anesthetic Gas Module) M1026A/B: 18 lbs (8.2 kg)* M1013A: 12 lbs (5.45 kg)* M1019A: 12 lbs (5.45 kg)*</p>	 <p>TTM-05-050</p>	N/A	
	 <p>RS6-420-PKG</p>	 <p>AHM-053 (with Philips supplied Tabletop Mounts)</p>	N/A	

*Weight provided to help approximate mounting capacity. For reference only.

Philips Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>Multi-Measurement Module (MMS) 1.4 lbs (0.65 kg)*</p>	N/A (Mounts to the back of the MP20/30)	N/A	N/A	N/A
 <p>Cable and Patient Leads</p>	 <p>CRDHLDR</p>	 <p>AHM-053</p>	Mounts below with primary device on AHM/SSM	N/A
		 <p>SSM-053</p>		
	 <p>101004-MON-BB</p>	N/A	Mounts on MRS rail with AHM/SSM	N/A

*Weight provided to help approximate mounting capacity. For reference only.