

Titan Bariatric Stretcher



Key Features:

- 700 lbs (318 kg) weight capacity
- Standard fifth wheel for improved mobility
- 21" (53 cm) low height
- Storage area for two oxygen E-cylinder tanks in base and patient belongings
- Radiolucent fowler section
- Dual brake pedal mechanisms
- Optional knee gatch for enhanced comfort
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces

The Amico Titan Bariatric Stretcher uses an innovative design and incorporates distinctive features aimed to help caregivers provide a more meaningful level of care to patients.

This product is specifically designed to meet the needs of customers by offering key characteristics such as a dual brake system, deep well storage area in the base and an ergonomic design to help improve efficiency of healthcare. With the heightened awareness regarding fall prevention, the Titan Bariatric Stretcher prides itself by offering the best-in-class low height stretcher on the market. This Amico product differentiates itself with a superior storage area in its base.

Standard Features:

- Fold-down siderails
- Four IV pole receptacles
- Roller bumpers
- Adjustable fowler section
- 8" (20 cm) wheel casters
- Fifth wheel steering
- Stationary foot end litter
- Hydraulic lift system capable of supporting up to 700 lbs (318 kg)
- Central brake locking system
- Ergonomic push handles
- Integrated foot push handle
- Side bumpers
- Head and foot-end mattress retainer

Model Number & Approval

Model Number	S-AM1-3003
Agency Approval	cULus

Overall Dimensions

Overall Length	84.5" (214 cm)
Overall Width	38" (96.5 cm)

Litter Positioning

Trend/Reverse Trend	17°
Maximum Fowler Section Inclination	80°
Maximum Knee Section Inclination	35°

Height Range

Lowest	21" (53 cm)
Highest	33" (84 cm)

Surface Dimensions

Surface Width	30" (76 cm)
Surface Length	75" (191 cm)
Surface Thickness	4" (10 cm)

Weight Capacity & Caster Dimension

Maximum Weight Limit	700 lbs (318 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- 3" (8 cm) ultra comfort mattress
- 4"/5" (10 cm/12.7 cm) high density foam mattress with or without fire barrier
- Foley bag holder
- Utility tray
- Permanent IV pole
- Overbed tables

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: Three years parts

Amico Beds Corporation 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 71 East Industry Court, Deer Park, NY 11729, USA
Toll Free Tel: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862