

Titan EYE/ENT/ORAL Care Surgery Stretcher



Key Features:

- 550 lbs (250 kg) weight capacity
- 21.25" (54 cm) low height
- Storage area for patient belongings
- Optional radiolucent fowler section
- Easily accessible brake pedals
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces
- Articulating headpiece

The Amico Titan EYE/ENT/ORAL Care Surgery Stretcher is designed with functionality and convenience in mind, allowing caregivers to provide a more meaningful level of care to patients.

This product offers key features such as an articulating headpiece that is a stainless steel design for enhanced motion, durability and sterilization; a deep well storage area for patient belongings and an overall ergonomic design to help improve efficiency of healthcare. To accommodate surgical patients with limited mobility, the Titan EYE/ENT/ORAL Care Surgery Stretcher is built on Amico's best-in-class, low height stretcher platform.

The Titan EYE/ENT/ORAL Care Surgery Stretcher can also serve as an interim transport stretcher for moving patients between departments

Standard Features:

- Four IV pole holders
- Roller bumpers
- Adjustable fowler section
- Fifth wheel steering
- Fold-down siderails
- 8" (20 cm) wheel casters
- Two side mounted pedals for patient surface height adjustment
- Hydraulic lift system capable of supporting up to 550 lbs (250 kg)
- Four-caster braking system
- Ergonomic push handles
- Side bumpers
- Head and foot-end mattress retainer
- Articulating headpiece

Model Number & Approval

Model Number	S-AM1-3004
Agency Approval	cULus

Overall Dimensions

Overall Length	86" (218 cm)
Overall Width	32.5" (82.6 cm)

Litter Positioning

Trend/Reverse Trend	±15.5°
Maximum Fowler Section Inclination	81.5°
Maximum Knee Section Inclination (optional)	24°

Height Range

Lowest	21.25" (54 cm)
Highest	33" (83.8 cm)

Headpiece

Height Adjustment	+4" (10 cm)
Side to Side Adjustment	±15°

Surface Dimensions

Surface Width	26" (66 cm)
Surface Length	83" (210 cm)
Surface Thickness	4" (10 cm)

Weight Capacity & Caster Dimension

Maximum Weight Limit	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- 3" (8 cm) ultra comfort mattress
- 4"/5" (10 cm/12.7 cm) high density foam mattress with or without fire barrier
- Foley bag holder
- Disposable waste basins
- Utility tray
- Permanent IV pole
- Overbed tables

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: Three years parts

Amico Beds Corporation 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 71 East Industry Court, Deer Park, NY 11729, USA
Toll Free Tel: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862