



Magna 22 is a conventional approach to magnification and inspection lighting. It provides 4 diopters of magnification and the most rugged construction of any "ring magnifier" in the industry. An optional snap-in Add-X lens increases magnification to 8 diopters. This light is ideal for dermatology and cosmetic applications.

## Features:

**Connected Load:** 120 V; 60Hz

**Fitted With:** 1 x fluorescent lamp FC8T9 22W; Neutral white

**Work Equipment:** Electronic ballast

**Mains Lead:** Approx. 6ft (1.8m)

**Luminaire body Material:** Plastic

**Tubular Sections Material:** Plastic

**Weight (NET):** Approx. 6.8lbs (3.1kg)

### Fastening:

- Table mount included for table mounted lights
- Rail mount included for rail mounted lights
- Wall mount included for wall mounted lights
- Floor stand included for floor stand light

**Technology:** Switchable

**Usage:** Switch

**Lamp Cover:** Polycarbonate; White opal

**Class Of Protection:** II

**Color Luminaire Body:** Light grey

**Color Of Tubular Sections:** Light grey

**Lens; Dimensions:**  $\varnothing$  4.7" (120mm)

**Lens; Diopters:** 4

**Add-On Lens:**  $\varnothing$  1.2" (30mm) / total 12 diopters

# Dimensions



## Model Examples

<b>L-MA22-AB-RH-A</b>	Magna 22, Articulating Arm, Basic, Rail Mount, Horizontal, Amico Rail
<b>L-MA22-AB-W</b>	Magna 22, Articulating Arm, Basic, Wall Mount
<b>L-MA22-AE-F</b>	Magna 22, Articulating Arm with Extension, Floor Stand
<b>L-MA22-AE-T</b>	Magna 22, Articulating Arm with Extension, Table Mount

**NOTE:** North American: 120 V; International: 230 V

## Ordering Guide

